



August 31, 2017 **BLOG**

## **XPure General Info**

The static SMB system [XPure-S](#) that represents the XPure patented valve system with a zero dead leg concept.

The carousel SMB system [XPure-C](#) featuring a turntable that carries the sorption containers going through the distinct process steps of a typical IX chromatography cycle.

XPure offers SMB technology for product recovery and purification, comprising of:

- Feasibility studies
- Developed simulation modeling software
- Full-scale design
- Start-up, validation, maintenance and (process) service

and

- Is not related to any manufacturer of resins
- Applies originally manufactured components or and fabricates that are consistent with clients material standards
- Optional elements
- Vessels and liquid distributors
- (Inter)connecting piping and auxiliary equipment

[Download brochure](#)